



Primary HPV Cervical Screening – A Laboratory Perspective

Dr Dave Nuttall

CervicalCheck Clinical Laboratory Advisor



An tSeirbhís Náisiúnta Scagthástála
National Screening Service


CervicalCheck
AN CLÁR NÁISIÚNTA SCAGTHÁSTÁLA CEIRBHEACS
THE NATIONAL CERVICAL SCREENING PROGRAMME



Learning Outcomes of the Presentation

To highlight and discuss the following:

- What is HPV, its carcinogenic role in cervical cancer and its detection as an application in Cervical Screening
- The importance of an evidence based approach
- The HPV test, including practical considerations and challenges about its implementation as a screening test
- The importance of Quality Assurance including validation and verification processes
- Limitations of the test including potential for false negatives
- Sample taker responsibilities and optimising laboratory screening services



Facts about Laboratory Cervical Screening

- Detect pre-cancerous conditions in a timely manner to prevent progression to invasive disease
- Reduce morbidity and ultimately mortality
- A HPV Test is a test of risk of cervical disease
 - Has a high sensitivity and negative predictive value
 - Poor specificity
- Therefore a cervical screening programme also requires a test of disease
 - High specificity
 - High positive predictive value
 - Cytology
 - Molecular Biomarkers
 - Computer Assisted Screening
 - Digital Artificial Intelligence Screening
- Current measure of outcome of a cervical screening programme is cervical histopathology



Following the Evidence

- 1842 – Domenico Antonio Rigoni-Stern – Statistical facts about cancers – IV Conference of Italian Scientists
- George Herbert Green, “The Unfortunate Experiment” and the Cartwright report
- 1976 – Harald Zur Hausen - HPV related to cervical cancer
- HPV testing initially introduced as a risk stratification tool
- HPV triage of low grade cytology
- Test of Cure following treatment
- Primary HPV test – numerous clinical trials - ARTISTIC, ATHENA
- Countries now implementing primary HPV screening include:
- Australia, Denmark, England, Holland, Ireland, Scotland, Wales



Health technology assessment of human papillomavirus testing as the primary screening method for prevention of cervical cancer

24 May 2017

Safer Better Care

Human Papilloma Virus - HPV

Over 100 sub-types, but only a relatively small number have been implicated in cervical disease




HPV mediated carcinogenesis

- HPV infection of epithelium
- Immune system deficiencies
- Expression of E6, E7 oncogenes
- Cellular alteration and transformation
- Dyskaryosis and malignancy



HPV Related Cervical Carcinogenesis



- 80% of infections are cleared within 2 years
- Lead time to developing cancer is 5-20 years
- Association with persistent infection



The challenges of HPV testing in cervical primary screening

HPV Not Detected Result

- Negative predictive value 99.9%
- Long-term protective effect is well established
- Safely increase the screening interval

HPV Detected Result

- Poor specificity of a positive HPV test
- Risk of over referral to colposcopy
- Increased cost and patient anxiety



Specificity
NPV

vs

Sensitivity
PPV



The HPV Test – Clinical Validation

- Currently over 200 available HPV assays
- CervicalCheck has a policy of rigorous validation of HPV tests
- Those that satisfy programme requirements are added to the Approved Assay List
- HPV detection threshold – approved tests have been clinically validated so that a positive or negative result can be trusted. There is no need to be cautious of results that are “close to the cut-off”
- This is why a negative result is has such a high negative predictive potential
- The test is very reliable, **IF** the CervicalCheck standards and processes are followed



Operational considerations of introducing primary HPV testing in the laboratory



Acknowledgments and the future!

- The dedication and professionalism of our laboratory providers
 - Quest Diagnostics
 - Coombe Women and Infants University Hospital
- National Cervical Screening Laboratory development
 - An exciting project and very positive for cervical screening in Ireland
 - Challenging – staffing in Ireland at critically low levels
 - Training infrastructure for Cytopathologists needed
 - Advanced clinical medical scientist practitioners role needed urgently



The HPV Test – General Information

- Platform currently in use – Roche Cobas 4800
- Tests for 14 high risk sub-types
- Partial genotyping results are not provided for screening samples
- Test for infection – therefore a test of **RISK** for cervical disease
- Negative HPV test means very, very small chance of abnormality
- Cytology is the test used as reflex to a HPV positive result. It is a test of **DISEASE**, specifically employed to qualify the **RISK**
- All high risk HPV positive samples have a reflex cytology test performed to check for the presence of cervical disease
- Automated and therefore has a faster TAT than a cytology primary test



Partial Genotyping



- Widely published that most (up to 70%) of HG cervical lesions attributable to HPV subtypes 16,18
- HPV testing demonstrates that HPV vaccination is performing, with declines in HPV 16 and 18 prevalence recorded
- This evidence^{1,2} has informed CervicalCheck primary HPV partial genotyping policy – reporting only presence or absence of high risk subtypes
- All positives referred for colposcopic assessment on second HPV positive
- Facilitates benchmarking with primary HPV based programmes e.g. UK
- Use of partial genotyping in colposcopy setting currently under review

1. Type-specific HPV in sexually active young females in England, to end 2018. Public Health England. PHE publications gateway number GW-1027; HPR 14(2). Available at: https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/858872/hpr0220_HPV_2018.pdf

2. Age-specific HPV prevalence among 116,052 women in Australia's renewed cervical screening program: A new tool for monitoring vaccine impact. Julia ML. Brotherton, David Hawkes, Farhana Sultana, Michael J. Malloy, Dorothy A. Machalek, Megan A. Smith, Suzanne M. Garland, Marion Saville. <https://doi.org/10.1016/j.vaccine.2018.11.075>



HPV Vaccinated Participants

- HPV vaccinated persons are increasingly tested as part of a cervical screening programme as they reach eligible screening age
- HPV vaccination results in fewer persistent infections, with subsequently fewer positive results
- Older HPV assays were not clinically validated in a vaccinated population and so further validation studies are required:
 - Performance evaluation of the assay concerned in other primary HPV based screening programmes where HPV vaccinated eligible participants are screened (Australia, U.K.) is important



Vaginal Vault Samples

- Traditionally, eligible post hysterectomy samples followed up with cytology
- Currently, no HPV tests are FDA approved for HPV testing on vaginal samples
- Non-FDA in-house analytical validation carried out – good concordance with cervical sample testing
- Clinical validation carried out on samples with known cytology outcomes
- Initial correlation results are promising
- Now looking to perform histology outcome correlation



The “HPV – indeterminate” result

This is a relatively rare event, however, there are a number of causes, including :

- Assay quality control failure
- Technical issues with the sample or equipment
- Sample contamination

In many instances, the test can be repeated on the sample and a result obtained. However, that is not always the case – hence the indeterminate result.



Sample Contamination

- All HPV assays can be affected by contamination which inhibits the test reaction, potentially resulting in a false negative result
- From the programme perspective, this mainly concerns blood and certain lubricants
- The tolerance for contamination with whole blood varies with assay, generally the newer assays have higher tolerance
- Important that the laboratory has a standardised, reproducible process that minimises uncertainty around the test



Role of the Doctor/Nurse

- Fundamental to the success of the screening programme – interfacing directly with the woman
- Poor sampling can lead to failure in detecting pre-malignant abnormalities - before the sample has even reached the lab!
- Sampling is somewhat less critical with primary HPV testing, but is still a major consideration for reflex cytology
- Fact - Abnormal results are relatively commonplace but cancer is relatively rare
It is vital, therefore, that the smear taker visualises the cervix when sampling - the lab can only give a result on the material presented in the vial
- **HELP THE LAB TO HELP YOU!**



Specimen Collection – Request Form

- Ensure that the current version of the form is used
- For electronic requests, ensure that the GP IT system provider has the current request order version
- Check that the patient details are correct – **WHILE THE PATIENT IS THERE!**
- Rejection by the laboratory because of incompletely completed request forms may result in the sample becoming date expired by the time a corrected form is received
- Ensure that all relevant clinical history is captured
- **CONSENT** – ensure that this is captured to latest requirements



Specimen collection – Date Expired Samples

- GLP requires that all lab investigations producing test results are controlled where possible, to ensure the most accurate results
- Laboratories therefore estimate and minimise uncertainty wherever and whenever possible
- The sample collection fluid is effective in preserving cells for 42 days, after which its properties begin to decline
- HPV is stable up to 6 months from taking the sample
- Because of the potential requirement for reflex cytology, a 42 day expiration limit is operated
- Samples in date expired vials will be rejected, ensure availability of in-date vials
- Expiry dates include the 42 day period, and so stock control in your areas is absolutely vital
- Do not delay vials to save on postage, could make the difference in none, some or even all samples being rejected and the woman having to attend for a repeat test



Specimen collection – taking the sample

- Bin spare labels and forms – minimise the chance of transcription/ID errors
- The sample life for a HPV test and subsequent cytology test (if required) is 42 days – this has implications for the sample taker
- Make sure that the vial is not date expired, nor likely to in the near future that might impact on the validity of the test
- Sample transport – ensure that the samples get to the laboratory in a timely manner, the 42 day clock is ticking as soon as the broom enters the sample collection fluid
- Batching samples together to save postage can have serious impact on the test validity



Common Reasons for Test Rejection by the Lab

- Patients outside programme eligible screening age
- Request form incomplete/incorrect version of the form
- Sample not labelled
- Consent to test not adequately recorded/no signature
- Sample vial date expired or within 42 days of expiry
- Sample leaked
- Broom left in vial



Tests Queried: October 1st to November 15th, 2020



Expired Vials	237	10.9%
Address Mismatch	8	0.4%
DOB mismatch	37	1.7%
DOS mismatch	69	3.2%
MCRN missing/incorrect	1551	71.1%
Blank Vial	84	3.8%
Name Mismatch	22	1.0%
Others	10	0.5%
Incorrect Forms	163	7.5%
TOTAL	2181	



Impact of Rejected/Queried Tests for Patient and Lab


- In 6 weeks, over 2000 samples were rejected or queried
- In total the lab would expect to receive about 36000 samples, so about 6% of samples are either queried or rejected
- Follow up of these samples has a major impact on available resources
- A significant number of women undergo a repeat sample
- Programme reputation potentially damaged with every incident of this nature
- The vast majority of these issues are AVOIDABLE!

The ThinPrep Liquid Based Cytology test






Cells.....!



HSIL – microbiopsy



LSIL – Koilocytes



Cytology and TZ Sampling

- Data will no longer be collected on Transformation Zone sampling as this will only be available on HPV positive results – up to 13% of the total workload
- Sample taker proficiency will be quality assured in other ways – approved training, sample taker register, updates and inadequacy rates



The Clinically Indicated Cervical Screening Sample

- Should not be taken - a clinical referral is more appropriate
- HPV primary screening test is **not** a diagnostic test (designed to detect HPV as causative to cervical pre cancer – can be negative where the HPV DNA is integrated into the host genome)
- Sample takers must act appropriately in managing symptomatic patients



The Private Test

A HPV or cytology test can be negative or positive for a number of reasons – not necessarily associated with cervical disease. Most labs are accredited in one form or another, however:

- CervicalCheck quality assures the labs it contracts with, which are then approved
- Labs are re-assessed if circumstances change and routinely within 2 - 3 years
- CervicalCheck cannot accept results from non-approved labs and when the programme is aware of the existence of such a result, further steps will be taken to assure the validity of that result
- CervicalCheck is currently reviewing the process for managing results from non-approved sources to ensure they are appropriately incorporated into the governance and assurance processes of the programme



Thanks for Listening!



An tSeirbhís Náisiúnta Scagthástála
National Screening Service


CervicalCheck
AN CLÁR NÁISIÚNTA SCAGTHÁSTÁLA CEIRBHEACS
THE NATIONAL CERVICAL SCREENING PROGRAMME