


## Guidance Note 2

### Using the peel-off vial number to help match vials

Each vial has a unique number. The number is printed twice on each vial label, one on the permanent section of the label and the second on the peel-off section of the label.


The peel-off vial number label adds an additional unique identifier to match the vial to the screening form during transit and at reception at the laboratory.

#### 1. Remove the peel-off unique number from the vial label



Note: The woman's name and DOB **must still be written** on the vial label. The vial number provides an additional printed identifier to help match vials to forms at the laboratory.

#### 2. Stick the unique vial number to the correct Cervical Cytology Form



Note: Please ensure you are using the latest revision of the Cervical Screening Form. You will find the most current version of the Cervical Screening Form on [www.hse.ie/cervicalcheck](http://www.hse.ie/cervicalcheck).

\*Please note, throughout this document, where we refer to 'women', we mean 'women, or people with a cervix'.